



2018


Conclusion Of Machine Learning Of Kaggle

Lorem ipsum dolor sit amet, consectetur adipiscing elit.

Succinct Introduce Competition

- Competition platform: Kaggle
- Sponsor Institution: RSNA (Radiological Society of North America)
- Data Furnish Institution: NIH (National Institutes of Health Clinical Center)
- Total Fund: \$30,000
 - First: \$12,000
 - Second: \$7,000
 - Third: \$4,000
 - Fourth to Tenth: Each \$1,000
- Total Joined Teams: 1445

Distribution Of Data



Final Rank

- The Competition Max Score is : 0.26
- Our Max Score is : 0.181
- The Max Rank Of us is : 160 th
- The Gross teams are : 1445
- The percentage of us is : Top 11%

Processing

1

Initial Dilemma

2

Frame Of Model

3

Deficient Of Model

4



ResNeXt-FPN Frame

ONE


Initial Dilemma

Initial Dilemma

- At the beginning of the competition, I have used the single model Mask-rcnn to detect the Lung Opacity area in the X-ray image. However, the final result of the single model was not excellent. There were numerous false positive samples in it. I have extrapolated the reason of the outcome. The conclusion of mine is because of the similarity of the X-ray images in different classes. With the observation of my eyes, I can not discern the **Lung Opacity** patient and the **Not normal but not Lung Opacity** patient.
- So, with the conclusion that I have proposed, I thought we must use the multiple model to suppress the false positive in the samples which will be used in the final detection model. In the theory, it can pronounced suppress the false positive.




Not Normal Not Lung Opacity




Lung Opacity

Initial Dilemma


- The second of dilemma at the beginning is the deficient of X-ray image . With the over exposure , It is so vague in one image that can not distinguish the margin of the lung .
- So , for surmount the plight, We have tried numerous image trans , like **Gamma trans** , **Liner trans** at all . At the end , the **histogram equalization** is the final trans that we have choised.
- And this is one of our data augmentation method .The others augmentation have displaied at the side of the screen.




Original Image
3d63ab2e-862a-4f11-82d3-b9f561e8ecec



Reorganized Image



Original After Trans Image
3d63ab2e-862a-4f11-82d3-b9f561e8ecec




Reorganized Image After
Trans

TWO

Frame Of Model

Frame Of Model


- For solve the plight that we have encountered at the first part ,i proposed the model which has two steps , **classification and detection** .
- There is three labels in data set , **Normal** , **Not Normal / No Lung Opacity** , **Lung Opacity** . Our terminal purpose is to exclude more negtive samples and guaranteeing all Lung Opacity samples are all in the samples which will be detected .
- Detection will be used to detect the samples which have the Lung Opacity label .



Frame Of Model (Classification)

- The task of classification has **three** steps .
- **ONE** : Use classification to discern the **Normal** and **Not Normal** samples .The samples **Not Normal** contain Not Normal / No Lung Opacity and Lung Opacity label samples .Training samples contain all labels in the data set (Training samples : 24650 , validate samples :1034) .For guaranteeing the Lung Opacity samples are all included in the next step , i modulated the **sensitive ratio to 0.99** .
- **TWO** : Use the classification to discern the samples which label is **Not Normal** . To distinguish the Not Normal / No Lung Opacity and Lung Opacity samples . Training samples (Training samples : 16368 , validate samples :892) only contain labels of Not Normal / No Lung Opacity and Lung Opacity .(**Sensitive ratio is 0.99**).
- **THREE** : For excluding more negtive samples in the label of **Lung Opacity** (There are numerous negtive samples in the label of Lung Opacity , because the sensitive of the above double classification are 0.99!), I design the final classification. Use the classification to discern the negtive samples which are in the label of **Lung Opacity** .Training samples (Training samples : 24650 , validate samples :1034) contain all labels in the data set .(**Sensitive ratio is 1.0**)

Classification work flow



Frame Of Model (Classification)

- All the classification are based on the model of **ResNet 50** or **ResNet 75**. The First Classification is based on the **ResNet 50** , because it is more simple than the next two model . The second and third classification is based on **ResNet 75** , its are more complex than the First model .
- For improving the accuracy of each steps , i use the **combination trick** . In each step , the classification unites with four different subclassification with different accuracy . And the final result combines the distinct four result .
- In each subclassification , the accuracy will effect the weight of this subclassification. The sum of weights of all subclassification is 1.0 .

```
TP is 606
TP_FN is 704
TP / TP_FN 0.8607954545454546
TN is 305
TN_FP is 330
TN / TN_FP 0.9242424242424242
ACC is 0.8924242424242425
```

Normal or Not Normal
Accuracy (First Step)

```
TP is 139
TP_FN is 198
TP / TP_FN 0.702020202020202
TN is 549
TN_FP is 694
TN / TN_FP 0.7910662824207493
ACC is 0.7449494949494949
```

Not Normal / No Lung
Opacity or Lung Opacity
Accuracy(Second Step)

Combination
Final
Classification.

```
TP is 496
TP_FN is 530
TP / TP_FN 0.9358490566037736
TN is 532
TN_FP is 637
TN / TN_FP 0.8351648351648352
ACC is 0.8849056603773585
```

Lung Opacity or Others
Accuracy(Third Step)

Frame Of Model (Classification)


- For **validating** the final result of the classification , i random select 1167 samples in the data set to attest the efficiency of the classification.
- The final result has displayed on the side of the screen . The **sensitive ratio is 0.996**. At this condition , the classification can remove **0.62** negative samples in this data set .

```
TP is 263
TP_FN is 264
TP / TP_FN 0.9962121212121212
TN is 680
TN_FP is 1096
TN / TN_FP 0.6204379562043796
ACC is 0.8068181818181818
All numbers of Lung opacity 679
```


Holistic classification ratio of sensitive
and Specificity

Frame Of Model (Detection)

- We have test two detection model , Mask-rcnn and Yolo3 .
- In this competition , the model of Mask-rcnn performed weaker than the Yolo3 . But in the theory , Mask-rcnn should perform more robust and accuracy than the model of Yolo3.
- The most essential question is that we can not train a fantastic weights of Mask-rcnn .Not matter what we have done , forllowing the training steps in the thesis or training model with all losses.
- The detection model of Yolo3 performs better than the Mask-rcnn in the validate data set .So , we choised Yolo3 finally .







Mask-rcnn



Yolo-V3

Frame Of Model (Result Display)



THREE

Deficient Of Model

Deficient Of Model

- **ONE** : The original format of the data set is DCM . it contains many useful information in this file , like gender , age , view position ,and so on . However , we have ommited those useful information and postulated all Lung Opacity images have the **same distribution** . So this many deter our to raise accuracy .
- **TWO** : We can not training a excellent weights for detection model .In the validate set , the detection model performed well , but in the test set , it performed worse than the validate set . it always ommits or detects error .
- **THREE** : I have googled some information with the Lung Opacity . In some way , doctor can not decree one patient who has Lung Opacity or not only through the image of X-ray . There are also another check , like Blood examination , CT scan at al.For this reason , it may thwart our to improve the accuracy of detection .

FOUR

ResNeXt-FPN Frame

ResNeXt-FPN Frame(Frame Work)







Figure 3. A building block illustrating the lateral connection and the top-down pathway, merged by addition.



ResNeXt-FPN Frame(Display)



ResNeXt-FPN Frame(Compare)

- The subclassifications (ResNet75) of Thrid are all training with the same training set as ResNeXt50-FPN. The validate set are also all same .
- As show on the side of the screen. The layers of subclassifications are all deeper than the ResNeXt50-FPN, but accuracy are all lower than ResNeXt50-FPN. And in the exam of mine , the sub path way of the layer of ResNeXt50 is only 4 (Because tensorflow can not build so prodigious Graph) , less 8 times than original thesis .

```
TP 134
TP_FN 146
TN 674
TN_FP 888
TP / TP_FN 0.9178082191780822
TN / TN_FP 0.759009009009009
ACC is 0.8356164383561644
```

ResNeXt50-FPN Accuracy

```
TP is 163
TP_FN is 221
TP / TP_FN 0.7375565610859729
TN is 724
TN_FP is 889
TN / TN_FP 0.8143982002249719
ACC is 0.7737556561085973
```

SubClassification of
Third Accuracy

```
TP is 213
TP_FN is 221
TP / TP_FN 0.9638009049773756
TN is 575
TN_FP is 889
TN / TN_FP 0.6467941507311586
ACC is 0.8031674208144797
```

SubClassification of
Third Accuracy

```
TP is 187
TP_FN is 221
TP / TP_FN 0.8461538461538461
TN is 672
TN_FP is 889
TN / TN_FP 0.7559055118110236
ACC is 0.8009049773755657
```

SubClassification of Third
Accuracy

```
TP is 205
TP_FN is 221
TP / TP_FN 0.9276018099547512
TN is 612
TN_FP is 889
TN / TN_FP 0.688413948256468
ACC is 0.8076923076923077
```

SubClassification of
Third Accuracy



2018

THANKS

Lorem ipsum dolor sit amet, consectetur adipisicing elit.

